

Fibronectin FN

DS-Pr-R-903-A/3

Fibronectin FN
Recombinant Human Fibronectin FN Protein

/ 50µg / GMP-TL903-0050
250µg / GMP-TL903-0250

E.coli
Cold-insoluble globulin
DNA Fibronectin P02751 His N
64.7kD
> 90% SDS-PAGE
< 0.1EU/µg
hiPS ED₅₀ 3.2µg/mL

CS1

rhFN CH-296 II
VLA-4 / VLA-5 Hannenberg
50%~75%

-20
12
PBS -20 6 -80

-20

24

50 13 1 3

400-010-5556

1. Williams CM, Engler AJ, Slone RD, Galante LL, Schwarzbauer JE (2008) Fibronectin expression modulates mammary epithelial cell proliferation during acinar differentiation. *Cancer Research* 68 (9): 3185–3192.
2. Mao Y, Schwarzbauer JE (2005) Fibronectin fibrillogenesis, a cell-mediated matrix assembly process. *Matrix Biology : Journal of the International Society for Matrix Biology* 24 (6): 389–399.
3. Erickson HP (2002) Stretching fibronectin. *Journal of Muscle Research and Cell Motility* 23 (5-6): 575–580.

2023 10 17

2025 02 26